

# LipoTrust™ EX Gene



## Contents

### LipoTrust™ EX Gene

#### Amount and Storage

Lyophilized powder containing 1 μ mol of cationic lipid suitable for 1 ml To be kept at 2-8 C

#### Description

**LipoTrust™ EX Gene** Reagent is a proprietary cationic liposome formulation that facilitates highly efficient delivery of **Plasmid DNA** to mammalian cells.

Lyophilized powder is stable so that customers can use this as long as customer keeps as instructed.

#### How to use

### Transfection Procedure

#### A-1. Preparation of liposome suspension

1. Remove the aluminum cap of vial by tool with extreme care and take off the rubber cap.
2. Add 1 mL of **nuclease free** water into the **LipoTrust™ EX Gene** vial and gently shake the vial until no powder is remained. Be kept for 10 minutes before use. Please store this suspension at 2-8 .

#### A-2. Transfection

Use this procedure to transfect Plasmid DNA into mammalian cells using **LipoTrust™ EX Gene** per 1 well in a 24-well format. For other formats, see the table in **Recommended Reagent Amounts and Volumes** for appropriate reagent amounts to add.

**Tip:** To reduce well-to-well variability when transfecting multiple replicates (triplicates), proportionally scale up the reagent volumes to form complexes, then aliquot an equal volume of complexes into each well.

#### STEP 1. Cell plating

One day before transfection, plate cells and culture them for 1 day in ordinary growth medium without antibiotics per well. Cells should be 60-80% confluent at the time of transfection. (ex. HeLa Cell:  $2-4 \times 10^4$  cells/well)

#### STEP 2. Preparation of Transfection Complex and Application to Cells

1. Add 0.8μg of **Plasmid DNA** in 50 μl of serum free medium previously prepared in a suitable microtube and mix well.
2. Add 2 μl of prepared **LipoTrust™ EX Gene** suspension (A-1) into 50μl of serum free medium previously prepared in a suitable microtube and mix well gently.
3. Mix liquid of **STEP 2-1** with that of **STEP 2-2** well and kept for 20 minutes to form **Transfection complex** at room temperature.  
\* Either liquid can be poured into the other microtube to form **Transfection complex** in this case.
4. Transfer all of this **complex** to well in which cells were plated (STEP 1) and mix well with gently rocking the plate back and forth.

#### STEP 3. Incubation

Incubate the cells at 37 for 18 to 48 hours.

Incubation time may be varied according to cells or Plasmid DNA so that previous test to examine appropriate conditions is recommended.

Recommended Reagent Amounts and Volumes per 1 well

STEP	1	2-A)		2-B)	
Plate/dish	Medium volume/well	Serum-free medium dilution for Plasmid DNA	Plasmid DNA <sup>*1</sup>	Serum-free medium dilution for LipoTrust™ EX Gene	LipoTrust™ EX Gene <sup>*1</sup>
96 well	100 μl	25 μl	0.2 μg	25 μl	0.5 μl
48 well	200 μl	50 μl	0.4 μg	50 μl	1 μl
24 well	500 μl	50 μl	0.8 μg	50 μl	2 μl
12 well	1ml	100 μl	1.6 μg	100 μl	4 μl
6 well	2.5 ml	250 μl	4 μg	250 μl	10 μl
10 cm	15 ml	1.5 ml	24 μg	1.5 ml	60 μl

<sup>\*1</sup> Optimum amount of **Oligonucleotides** and **LipoTrust™ EX Gene** may be varied according to kind of cells or conditions. Variation from 1/4 to 4 times of **LipoTrust™ EX Gene** is recommended to obtain optimum transfection efficiency.

#### TIP for Transfection

1. Although good transfection efficiency is observed with **LipoTrust™ EX Gene** under serum containing medium (Ex:FEB 10%), use serum free medium at the time of preparation of **Transfection complex**.
2. Use healthy cells
3. Use medium free of antibiotics
4. Use high quality of oligonucleotide
5. Optimum amount of **Plasmid DNA** and **LipoTrust™ EX Gene** may be varied according to kind of cells or conditions. Variation from 1/4 to 4 times of **LipoTrust™ EX Gene** is recommended to obtain optimum transfection efficiency.

**This product is distributed for laboratory research use only. Caution: Not for diagnostic, food, cosmetic or household use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. Never inject human.**

#### **Terms and Condition of Use**

Hokkaido System Science Co., Ltd. (HSS) grants to you, its customer, a non-exclusive, non-transferable license to access the **LipoTrust™ EX Gene** transfection protocol applicable to HSS products, as set forth above (the "Protocol"). You acknowledge that the Protocol is protected by copyright. All rights not specifically granted to you are expressly reserved to HSS.

You may use Protocol for personal, educational, or scientific research or professional use, but in no case for a fee. You may not remove, obscure or modify any copyright or proprietary notices, author attribution or any disclaimer contained in the Protocol.

HSS PROVIDES THE PROTOCOL "AS IS", WITHOUT WARRANTIES OF ANY KIND, EITHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF TITLE, OR IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT; YOUR USE OF THE PROTOCOL IS AT YOUR OWN RISK; NEITHER HSS NOR ITS AFFILIATES, EMPLOYEES, OFFICERS OR DIRECTORS SHALL BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, CONSEQUENTIAL OR PUNITIVE DAMAGES ARISING OUT OF THIS LICENSE OR ANY USE OF THE PROTOCOL.


Some jurisdictions do not allow the exclusion or limitation of liability for certain damages, so the above limitation may not apply to you to the extent prohibited by such local laws; if so, then HSS's liability for damages hereunder shall not exceed an equal to the amounts paid by you hereunder, or one hundred dollars (\$100.00), whichever is less.

#### **Confidentiality**

You acknowledge the Protocol is HSS's confidential information and constitutes HSS's proprietary and valuable trade secrets. You agree to take reasonable care to protect the Protocol from disclosure to third parties, except as may be expressly authorized by HSS.

**For any inquiry about the product, please contact to**

Hokkaido System Science Co., Ltd.  
2-1, Shinkawa Nishi 2-1, Kitaku  
Sapporo, Hokkaido, Japan 001-0932  
TEL: 81-11-768-5901 FAX: 81-11-768-5951  
E-mail: bio@hssnet.co.jp  
URL: <http://www.hssnet.co.jp>



- Inner box can be used as a vial rack.

- Product appearance is subject to changes without notice.