

# LaboPass™ cDNA Synthesis Kit

Cat. No. CMRTK001 50 reactions

Cat. No. CMRTK002 100 reactions

Lot No.

**Note: For laboratory use only**

## Description

- The LaboPass™ cDNA synthesis kit is optimized for efficient and reproducible synthesizing first-strand cDNA from total RNA or poly(A) RNA. This kit provides the all components necessary for synthesis of cDNA including reverse transcriptase, dNTPs, reaction buffer, RNase inhibitor, primer and nuclease-free water.
- Included M-MuLV Reverse Transcriptase in this kit is a recombinant form of the reverse transcriptase from the Moloney Murine Leukemia Virus (M-MuLV) which possesses enhanced cDNA synthesis activity and reduced RNase H activity. Reduction of RNase H activity enables higher yield of full-length cDNA transcripts and increased thermostability.

## Contents

	50 rxn	100 rxn
• Reverse Transcriptase (200 U/μl)	50 μl	100 μl
• 5X RT Buffer	200 μl	400 μl
• dNTP (each 10 mM)	50 μl	100 μl
• RNase inhibitor (40 U/μl)	50 μl	100 μl
• Oligo (dT) <sub>18</sub> (100 μM)	50 μl	100 μl
• Random hexamer (0.2 μg/μl)	50 μl	100 μl
• Nuclease-free water	1 ml	1.5 ml

Store at -20°C

## Applications

- cDNA synthesis for RT-PCR and RT-qPCR
- cDNA synthesis for Cloning

## Storage and Stability

LaboPass™ cDNA synthesis kit is stable for 1 year when stored at -20°C.

## Quality Control

Each lot of Reverse Transcriptase, reaction buffer, RNase inhibitor and dNTPs is tested for contamination such as nuclease.

## General cDNA synthesis protocol

1. Add the following components into a sterile, nuclease-free tube for each 20 μl reaction :

- |                               |             |  |
|-------------------------------|-------------|--|
| • RNA                         | 1-8 μl      | 0.1 ng ~ 5 μg total RNA or 1 ng ~ 500 ng poly(A) RNA |
| • Oligo (dT) or Random primer | 1 μl        |  |
| • dNTP (each 10 mM)           | 1 μl        |  |
| • Nuclease-free water         | Up to 10 μl |  |

2. Heat for 5 min at 65°C, quick chill on ice and briefly spin-down.

3. Add the following components to the above mixture and gently mix :

- |                                |             |
|--------------------------------|-------------|
| • 5X RT reaction buffer        | 4 μl        |
| • RNase inhibitor (40 U/μl)    | 1 μl        |
| • M-MuLV Reverse Transcriptase | 1 μl        |
| • Nuclease-free water          | Up to 10 μl |

4. Incubate at 42°C, for 30 ~ 60 min.

5. Inactivate enzyme by heating at 70°C for 15 min.

6. Store at -20°C until ready for use.