

LaboPass™ FFPE RNA Isolation Kit

Product description

LaboPass™ FFPE RNA Isolation kit provides the simple and fast method for purification of total RNA from various biological samples, such as FFPE Tissue section, cultured cells and tissue. Highly pure RNA can be isolated in less than 30 minutes without any use of hazardous organic solvents and alcohol precipitation. Purified RNA can be directly used in various downstream applications without any further manipulations.

Kit contents

Contents*			Storage conditions
Cat. No.	CMFR0050	CMFR0100	
No. of preparation	50	100	Room Temperature
Spin column CD mini w/tube	50 ea	100 ea	
Buffer DRB	5 mL	10 mL	
Buffer DSB	1 mL	1 mL	
Buffer RAL**	40 mL	70 mL	
Buffer RW	40 mL	80 mL	
Buffer RWA (conc.)***	15 mL **	25 mL**	
Nuclease Free Water	15 mL	15 mL	
DNase I (Lyophilized)****	220 U	440 U	

* This kit is delivered under ambient condition. When being used immediately on arrival, all of components can be stored room temperature. But if the kit is going to be stocked for a long time, lyophilized DNase I should be stored at 2 ~ 8°C for optimal conservation. Long exposure to heat source can deteriorate the performance of kit significantly.

** During shipment or storage under cold ambient condition, a precipitate can be formed in buffer RAL. Heat the bottle at 20°C ~40°C to dissolve completely before use.

*** Buffer RWA is provided as concentrate. Ethanol must be added before first use as the indication on the bottle labels.

**** After reconstitution, DNase I solution should be stored at -20°C for optimal conservation of activity.

Product use limitations

LaboPass™ FFPE RNA Isolation Kit is intended for research uses only. This kit is not intended for diagnosis or treatment for human. All due care and attention should be exercised in the handling of the products.

Safety information

LaboPass™ FFPE RNA Isolation Kit contains irritants which are harmful when in contact with skin or eyes or inhaled or swallowed. Care should be taken when handling this product. Always wear gloves and eye protection and follow standard safety precautions.

Buffer RAL and buffer RW contain chaotropes which can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

Preventing RNase Contamination

RNase can be introduced accidentally during RNA preparation. Various factors such as reagents, air, dust, and human hands or skin, can be the sources of RNase contamination. Always wear disposable gloves and use sterile, disposable plastic wares. Use the pipette reserved for RNA work. Maintain a separate area for RNA work. Carefully clean the surfaces.

Required Materials Not Provided

Reagents	<ul style="list-style-type: none">- β-mercaptoethanol, ACS grade or better- 70% ethanol, ACS grade or better- Absolute ethanol, ACS grade or better- Hexadecane (deparaffination solution)
Disposable materials	<ul style="list-style-type: none">- RNase-free pipette tips- Sterile 1.5 mL microcentrifuge tubes
Equipments	<ul style="list-style-type: none">- Equipment for homogenizing sample- Devices for liquid handling- Microcentrifuge- Suitable protector (ex; lab coat, disposable gloves, goggles, etc)

Preparation of ethanol-added buffer

Buffer RWA is provided as concentrate. Absolute ethanol (ACS grade or better) should be added before first use as below.

Buffer name	Cat. No.	Volume of contents	Ethanol to be added	The final volume
RWA	CMFR0050	15 mL	60 mL	75 mL
	CMFR0100	25 mL	100 mL	125 mL

Reconstitution of DNase I

Before experiment, lyophilized DNase I should be reconstituted in buffer DSB. To obtain a working solution (2 U/uL), add 110 uL (#CMFR0050) or 220 uL (#CMFR0100) of buffer DSB (provided) to the tube containing a lyophilized DNase I. Do NOT vortex while dissolving. Store the reconstituted DNase I solution at -20°C.

DNA-free RNA

Most of DNA in the starting sample is removed during the preparation procedure. But, if DNA should be further eliminated from the preparation, it can be removed by treating the eluate with DNase I. Refer to the 'Appendix I'.

Protocols for the preparation of total RNA

A. Procedure for FFPE Tissue section

To check before start,

- This protocol is suitable for fresh, frozen, or reagent-stabilized tissue sample.
 - Prepare 1% β -mercaptoethanol(BME) in buffer RAL freshly. Otherwise, you can add 1% BME into buffer RAL on every single preparation.
 - Prepare 70% ethanol and absolute ethanol.
 - If the treatment of DNase I is not necessary, carefully read the procedure before start.
 - Prepare 'DNase I reaction mixture' as below; Make this mixture as just before step 12 as possible.
- ① Mix DNase I solution with buffer DRB for a preparation as below table.

DNase I stock conc.	DNase I sol. to be added per prep.	Buffer DRB to be added per prep.	Final volume of mixture
2 U/uL*	2 uL	68 uL	70 uL

* Concentration of provided DNase I

- ② Mix gently by pipetting without vortex.
- ③ Keep the mixture on ice until use.

- 1. Using a scalpel, trim excess paraffin off the sample block.**
- 2. Cut section 5~20 um thick.**
- 3. Immediately place the section in a 1.5 ml microcentrifuge tube(not provide), and close the lid.**
- 4. Add 320 uL Hexadecane (deparaffination solution, not provide), vortex for 10sec, centrifuge briefly to bring the sample to the bottom of the tube.**
- 5. Incubate at 56°C for 3 min, then remove supernatant. Be careful not to remove tissue sample.**
- 6. Homogenization : Add 350 uL buffer RAL(including 1% BME), then vortex for 30 sec.**
 - *It is essential to homogenize thoroughly the tissue in buffer RAL and to lyse completely the sample.*
 - *When using the fiber-rich tissue samples like heart and muscle, it is strongly recommended to apply lesser sample than 10 mg. The larger sample may not be lysed completely, and this will bring up poor result.*
 - *Use a double volume (700 uL) of buffer RAL for the larger sample than 20 mg, but it should not exceed 30 mg per preparation.*
- 7. Centrifuge at 10,000 xg for 2 mins and transfer carefully the supernatant to a new 1.5 mL microcentrifuge tube (not provided).**
 - *This step helps remove the debris of the homogenate and not clog a column membrane by the debris.*
- 8. Add 250 uL of 100% ethanol to the lysate and mix well by pipetting.**
 - *Do not centrifuge the ethanol mixture.*
- 9. Transfer 600 uL of the mixture into a spin column CD and centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.**
 - *If the mixture remains, repeat step 4 with the remaining mixture.*
 - *Make sure that any residual mixture does not remain in the column after centrifugation. If a residual mixture has*

remained, centrifuge again at full speed until all the solution has passed through completely.

10. If DNase I-treatment is not necessary, go to step 13.

11. Apply 350 uL of buffer RW to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

12. Apply 70 uL of DNase I reaction mixture on to the center of the spin column membrane and incubate for 10 mins at room temperature.

➤ *Refer to the section of 'To check before start' for the preparation of DNase I reaction mixture*

13. Apply 350 uL of buffer RW to the spin column and let it stand for 1 min. Centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

➤ *If DNase I-treatment is skipped, apply 600 uL of buffer RW to the spin column instead of 350 uL.*

14. Apply 600 uL of buffer RWA to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

15. Apply 600 uL of absolute ethanol to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

16. Centrifuge at full speed (>12,000 xg) for 1 min for removing the residual wash buffer and transfer the spin column into a new 1.5 mL microcentrifuge tube.

➤ *Residual ethanol may interfere with downstream applications. Care must be taken not to be contaminated by the carryover of ethanol.*

17. Apply 50 uL of Nuclease Free Water to the center of spin column membrane and let it stand for 1 min.

➤ *Elution volume can be adjusted according to the purpose of experiment.*

➤ *Using lesser amount than 50 uL will decrease the total RNA yield but increase the concentration of RNA. Note that at least 30 uL of eluent must be applied because the lesser eluent will not soak the membrane entirely, followed by poor or unpredictable result.*

18. Centrifuge at full speed for 1 min for eluting the RNA.

➤ *Purified RNA can be temporarily stored at 4°C for an immediate analysis, but it is strongly recommended to store at -70°C for long-term storage.*

B. Procedure for animal cultured cells

To check before start,

- Prepare 70% ethanol and absolute ethanol.
 - Reconstitute the lyophilized DNase I as described in previous page.
 - Prepare 'DNase I reaction mixture' as below; *Make this mixture as just before step 7 as possible.*
- ① Mix DNase I solution with buffer DRB for a preparation as below table.

DNase I stock conc.	DNase I sol. to be added per prep.	Buffer DRB to be added per prep.	Final volume of mixture
2 U/uL*	2 uL	68 uL	70 uL

* The concentration of the provided DNase I

- ② Mix gently by pipetting without vortex.
 - ③ Keep the mixture on ice until use.
- If the treatment of DNase I is not necessary, carefully read the procedure before start.

1. Harvest cells in a tube as below.

➤ *Cells grown in monolayer*

Harvest up to 5×10^6 cells carefully using scraper, pellet cells by centrifugation at low speed (below 800 xg) for 5 minutes, and then discard the culture medium.

➤ *Cells grown in suspension*

Pellet up to 5×10^6 cells by centrifugation at low speed (below 800 xg) for 5 minutes and discard the culture medium.

➤ *Do not wash the cells before lysing with buffer RAL as this may cause mRNA degradation.*

2. Apply 350 uL of buffer RAL to the tube and make the mixture homogenate completely by pipetting or micro-homogenizer.

➤ *It is essential for good result to completely make the lysate homogenate. Incomplete homogenization can lead to a significantly reduced yield.*

➤ *If the number of cells exceeds 5×10^6 , add 700 uL of buffer RAL. Do not apply more than 1×10^7 cells for 1 preparation. An insufficient lysis due to large sample may result in poor yields.*

3. Add 1 volume of 70% ethanol to the lysate and mix thoroughly by pipetting.

➤ *Do not centrifuge the ethanol mixture.*

4. Transfer 700 uL of the mixture into a spin column CD and centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

➤ *If the mixture remains, repeat step 4 with the remaining mixture.*

5. If DNase I-treatment is not necessary, go to step 8.

6. Apply 350 uL of buffer RW to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

7. Apply 70 uL of DNase I reaction mixture on to the center of the spin column membrane and incubate for 10 mins at room temperature.

➤ *Refer to the section of 'To check before start' for the preparation of DNase I reaction mixture*

8. Apply 350 uL of buffer RW to the spin column and let it stand for 1 min. Centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

➤ *If DNase I-treatment is skipped, apply 600 uL of buffer RW to the spin column instead of 350 uL.*

9. Apply 600 uL of buffer RWA to the spin column, centrifuge at 8,000 xg for 30 secs, empty the

collecting tube and re-insert the spin column back to the tube.

10. Apply 600 uL of absolute ethanol to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

11. Centrifuge at full speed (>12,000 xg) for 1 min for removing the residual wash buffer and transfer the spin column into a new 1.5 mL microcentrifuge tube.

➤ *Residual ethanol may interfere with downstream applications. Care must be taken not to be contaminated by the carryover of ethanol.*

12. Apply 50 uL of Nuclease Free Water to the center of spin column membrane and let it stand for 1 min.

➤ *Elution volume can be adjusted according to the purpose of experiment.*

➤ *Using lesser amount than 50 uL will decrease the total RNA yield but increase the concentration of RNA. Note that at least 30 uL of eluent must be applied because the lesser eluent will not soak the membrane entirely, followed by poor or unpredictable result.*

13. Centrifuge at full speed (>12,000 xg) for 1 min for eluting the RNA.

➤ *Purified RNA can be temporarily stored at 4°C for an immediate analysis, but it is strongly recommended to store below -70°C for long-term storage.*

C. Procedure for animal tissues

To check before start,

- This protocol is suitable for fresh, frozen, or reagent-stabilized tissue sample.
 - Prepare 1% β -mercaptoethanol(BME) in buffer RAL freshly. Otherwise, you can add 1% BME into buffer RAL on every single preparation.
 - Prepare 70% ethanol and absolute ethanol.
 - If the treatment of DNase I is not necessary, carefully read the procedure before start.
 - Prepare 'DNase I reaction mixture' as below; Make this mixture as just before step 7 as possible.
- ④ Mix DNase I solution with buffer DRB for a preparation as below table.

DNase I stock conc.	DNase I sol. to be added per prep.	Buffer DRB to be added per prep.	Final volume of mixture
2 U/uL*	2 uL	68 uL	70 uL

* Concentration of provided DNase I

- ⑤ Mix gently by pipetting without vortex.
- ⑥ Keep the mixture on ice until use.

1. Homogenize up to 20 mg of animal tissue as one of the methods described below.

- *It is essential to homogenize thoroughly the tissue in buffer RAL and to lyse completely the sample.*
 - *When using the fiber-rich tissue samples like heart and muscle, it is strongly recommended to apply lesser sample than 10 mg. The larger sample may not be lysed completely, and this will bring up poor result.*
 - *Use a double volume (700 uL) of buffer RAL for the larger sample than 20 mg, but it should not exceed 30 mg per preparation.*
- A. Grind the tissue sample to a fine powder with liquid nitrogen in a pre-chilled mortar and pestle. Put up to 20 mg of the powdered tissue into 1.5 mL microcentrifuge tube. Add 350 ul of buffer RAL (including 1% BME) and pulse-vortex for 30 secs.**
- B. Homogenize up to 20 mg of the tissue sample in 350 ul of buffer RAL (including 1% BME) using a hand-held or a rotor-stator homogenizer.**
- C. Pulverize up to 20 mg of the tissue sample in 2.0 mL tube using bead-beater. Add 350 ul of buffer RAL (including 1% BME) and pulse-vortex for 30 secs.**
Follow the instruction manuals for usage of the bead-beating instruments.

2. Centrifuge at 10,000 xg for 2 mins and transfer carefully the supernatant to a new 1.5 mL microcentrifuge tube (not provided).

- *This step helps remove the debris of the homogenate and not clog a column membrane by the debris.*

3. Add 1 volume of 70% ethanol to the lysate and mix well by pipetting.

- *Do not centrifuge the ethanol mixture.*

4. Transfer 700 uL of the mixture into a spin column CD and centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.

- *If the mixture remains, repeat step 4 with the remaining mixture.*
- *Make sure that any residual mixture does not remain in the column after centrifugation. If a residual mixture has remained, centrifuge again at full speed until all the solution has passed though completely.*

5. If DNase I-treatment is not necessary, go to step 8.

- 6. Apply 350 uL of buffer RW to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.**
- 7. Apply 70 uL of DNase I reaction mixture on to the center of the spin column membrane and incubate for 10 mins at room temperature.**
 - *Refer to the section of 'To check before start' for the preparation of DNase I reaction mixture*
- 8. Apply 350 uL of buffer RW to the spin column and let it stand for 1 min. Centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.**
 - *If DNase I-treatment is skipped, apply 600 uL of buffer RW to the spin column instead of 350 uL.*
- 9. Apply 600 uL of buffer RWA to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.**
- 10. Apply 600 uL of absolute ethanol to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.**
- 11. Centrifuge at full speed (>12,000 xg) for 1 min for removing the residual wash buffer and transfer the spin column into a new 1.5 mL microcentrifuge tube.**
 - *Residual ethanol may interfere with downstream applications. Care must be taken not to be contaminated by the carryover of ethanol.*
- 12. Apply 50 uL of Nuclease Free Water to the center of spin column membrane and let it stand for 1 min.**
 - *Elution volume can be adjusted according to the purpose of experiment.*
 - *Using lesser amount than 50 uL will decrease the total RNA yield but increase the concentration of RNA. Note that at least 30 uL of eluent must be applied because the lesser eluent will not soak the membrane entirely, followed by poor or unpredictable result.*
- 13. Centrifuge at full speed for 1 min for eluting the RNA.**
 - *Purified RNA can be temporarily stored at 4°C for an immediate analysis, but it is strongly recommended to store at -70°C for long-term storage.*

Trouble shooting guide

Facts	Possible causes	Recommendations
Low recovery of RNA	<i>Sample not homogenized completely</i>	Incomplete homogenization can lead to decrease in yield of total RNA and this may be attributed to several reasons; - Insufficient mixing with buffer RAL - Too much mass in the starting sample - Poor pulverization of sample Ensure the complete homogenization of the sample with buffer RAL.
	<i>Too much starting materials</i>	Too much starting materials will bring about inefficient lysis, followed by poor RNA yields. Reduce the amount of starting material as described on procedure.
	<i>Poor quality of starting material</i>	Use a freshly harvested sample if possible. The harvested sample should be handled under low temperature before addition of buffer RAL.
	<i>Ethanol was not added to the wash buffer</i>	Buffer RWA is provided as concentrate. Ethanol must be added to this buffer before first use. If not, the result will be significantly poor.
	<i>Culture media not completely removed</i>	Culture media affect the lysis and the binding efficiency. Discard the culture media as completely as possible when harvest.
Column clogged	<i>Sample not homogenized completely</i>	Refer to the same item at 'Low recovery of RNA' section.
	<i>Too much starting materials</i>	Refer to the same item at 'Low recovery of RNA' section.
RNA degraded	<i>Inappropriate handling of starting materials</i>	Starting sample should be quickly treated under low temperature. Higher temperature or retarded processing would be a cause of degradation.
	<i>Poor quality of starting material</i>	Refer to the same item at 'Low recovery of RNA' section.
	<i>RNase contamination</i>	RNase can be introduced accidentally into a preparation at any steps. Always wear disposable gloves and use RNase-free plasticwares. Do not use shared equipment if possible.
DNA contamination	<i>Incorrect treatment of DNase I reaction mixture</i>	DNase I reaction mixture should be pipetted onto the center of the spin column membrane for proper enzymatic reactions.
Enzymatic reaction is not performed well with purified RNA	<i>Salt carryover in eluate</i>	Ensure that washing steps are carried out just in accordance with the protocols. Additional washing steps with absolute ethanol may help remove salts from the membrane.
	<i>Residual ethanol in eluate</i>	The spin column membrane should be dried completely before eluting. Perform additional centrifugation to dry the membrane, if needed. Do NOT incubate the column at high temperature.

Appendix I.

Procedure for DNase I treatment in RNA eluate

To check before start,

- Prepare 0.25M EDTA and 75°C water bath or dry bath.

1. Prepare the reaction mixture freshly in a 1.5 mL tube as below;

- 50 uL RNA eluate
- 5 uL buffer DRB
- 2U DNase I solution

2. Incubate the mixture for 10 mins at room temperature.

3. Proceed to next step for direct use or proceed to 'Appendix II' for clean-up the reaction mixture

- *It is recommended to proceed to clean-up procedure at 'Appendix II', because the inactivation process of DNase I at step 5 can deteriorate the quality of RNA.*

4. Add 1 uL of 0.25 M EDTA to the tube.

5. Inactivate DNase I by incubating at 75°C for 10 mins.

Appendix II.

Procedure for RNA clean-up

To check before start,

1. Adjust the sample volume to 100 uL with Nuclease Free Water.
2. Add 350 uL of buffer RAL and mix thoroughly.
3. Add 250 uL of absolute ethanol to the tube and mix well by pipetting.
4. Transfer all of the mixture into a mini spin column and centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.
5. Apply 600 uL of buffer RWA to the spin column, centrifuge at 8,000 xg for 30 secs, empty the collecting tube and re-insert the spin column back to the tube.
6. Repeat the step 5 once.
7. Centrifuge at full speed (>12,000 xg) for 1 min for removing the residual wash buffer and transfer the spin column into a new 1.5 mL microcentrifuge tube.
 - Residual ethanol may interfere with downstream applications. Care must be taken not to be contaminated by the carryover of ethanol.
8. Apply 50 uL of Nuclease Free Water to the center of spin column membrane and let it stand for 1 min.
 - Elution volume can be adjusted according to the purpose of experiment. Using lesser amount than 50 uL will decrease the total RNA yield but increase the concentration of RNA. Note that at least 30 uL of eluent must be applied because the lesser eluent will not soak the membrane entirely, followed by poor or unpredictable result.
9. Centrifuge at full speed for 1 min for eluting the RNA.
 - Purified RNA can be temporarily stored at 4°C for an immediate analysis, but it is strongly recommended to store at -70°C for long-term storage.