

ATTO-TEC 蛍光色素 一覧 (全33種類)

No.	蛍光色素名	最大吸収波長 (nm)	最大蛍光波長 (nm)	対応クエンチャー (nm)				
1	ATTO 390	390	476	Dabcyl 400-530				
2	ATTO 425	439	485					
3	ATTO 465	453	506					
4	ATTO 488	500	520					
5	ATTO 495	498	526					
6	ATTO 514	511	532		TAMRA 520-560	BHQ-1 480-580		
7	ATTO 520	517	538					
8	ATTO 532	532	552					
9	ATTO Rho6G	533	557					
10	ATTO 542	542	562					
11	ATTO 550	554	576					
12	ATTO Rho3B	566	589					
13	ATTO 565	564	590					
14	ATTO Rho11	572	595					
15	ATTO Rho12	576	601					
16	ATTO Thio12	582	607					
17	ATTO Rho101	587	609					
18	ATTO 590	593	622				BHQ-2 560-670	
19	ATTO 594	603	626					
20	ATTO Rho13	603	627					
21	ATTO 610	616	633					
22	ATTO 620	620	642					
23	ATTO Rho14	626	646					
24	ATTO 633	630	651					
25	ATTO 647N	646	664					
26	ATTO 647	647	667					
27	ATTO 655	663	680					
28	ATTO 665	662	680					
29	ATTO Oxa12	662	681					
30	ATTO 680	681	698					
31	ATTO 700	700	716					
32	ATTO 725	728	751					
33	ATTO 740	743	763				BHQ-3 620-730	